
Product Name: PCYT1A Rabbit Monoclonal Antibody**Catalog #: AMRe84553**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC,FC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.71mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
Molecular Weight	42 kDa

Antigen Information


Gene Name	PCYT1A
Alternative Names	CCT A; CCT alph; CTPCT; PCYT1; Pcyt1a;;PCYT1A
Gene ID	
SwissProt ID	P49585
Immunogen	A synthesized peptide derived from human PCYT1A

Background

Catalyzes the key rate-limiting step in the CDP-choline pathway for phosphatidylcholine biosynthesis.

Research Area

Image Data



Western blot analysis of PCYT1A expression in A549 cell lysate.

All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

