
Product Name: Acetyl-Histone H3(Lys27) Rabbit Monoclonal Antibody**Catalog #: AMRe84555**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,ICC,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Acetylated
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,IP 1:20-1:50
Molecular Weight	Calculated MW: 15 kDa ; Observed MW: 17 kDa

Antigen Information

Gene Name	Histone H3(acetylK27)
Alternative Names	Histone H3.1, Histone H3, HIST1H3A;;Acetyl-Histone H3 (K28)
Gene ID	
SwissProt ID	P68431/P84243/Q16695/Q6NXT2/Q71DI3
Immunogen	A synthesized peptide derived from human Histone H3.1 around the acetylation site of K28


Background

Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA


repair, DNA replication and chromosomal stability. DNA accessibility is regulated via a complex set of post-translational modifications of histones, also called histone code, and nucleosome remodeling.

Research Area


Image Data



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of Histone H3 (acetyl K27) expression in HeLa cell treated with TSA cell lysate.