

Product Name: SynapsinII Rabbit Monoclonal Antibody
Catalog #: AMRe84557



Summary

Production Name	SynapsinII Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	SynapsinII
Alternative Names	SYN2; Synapsin 2; Synapsin2; SynapsinII; SYNII; SYNIIa; SYNIIb;;Synapsin 2
Gene ID	
SwissProt ID	Q92777. A synthesized peptide derived from human Synapsin 2

Application

Dilution Ratio	WB:1:1000-1:2000
Molecular Weight	63 kDa

Background

Neuronal phosphoprotein that coats synaptic vesicles, binds to the cytoskeleton, and is believed to function in the


Product Name: SynapsinII Rabbit Monoclonal Antibody
Catalog #: AMRe84557



regulation of neurotransmitter release. May play a role in noradrenaline secretion by sympathetic neurons.

Research Area

Image Data



Western blot analysis of Synapsin II expression in Neuro-2a cell lysate.

Note

For research use only.