

Product Name: BMP15 Rabbit Monoclonal Antibody**Catalog #: AMRe84567**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,FC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,FC 1:20-1:100
Molecular Weight	Calculated MW: 45 kDa ; Observed MW: 50 kDa

Antigen Information


Gene Name	BMP15
Alternative Names	BMP15; GDF9B; Growth/differentiation factor 9B; ODG2; POF4;;BMP15
Gene ID	
SwissProt ID	O95972
Immunogen	A synthesized peptide derived from human BMP15

Background

May be involved in follicular development. Oocyte-specific growth/differentiation factor that stimulates folliculogenesis and granulosa cell (GC) growth.

Research Area

Image Data



Western blot analysis of BMP15 expression in HeLa cell lysate.