

Product Name: DP1 Rabbit Monoclonal Antibody**Catalog #: AMRe84581**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,FC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,FC 1:20-1:100
Molecular Weight	Calculated MW: 45 kDa ; Observed MW: 50 kDa

Antigen Information

Gene Name	DP1
Alternative Names	DP1; DRTF1; TFDP1;;DP1
Gene ID	
SwissProt ID	Q14186
Immunogen	A synthesized peptide derived from human DP1


Background

Can stimulate E2F-dependent transcription. Binds DNA cooperatively with E2F family members through the E2 recognition site, 5'-TTTC[CG]CGC-3', found in the promoter region of a number of genes whose products are involved in cell cycle regulation or

in DNA replication

Research Area

Image Data



Western blot analysis of DP1 expression in Raji cell lysate.