

Product Name: HMGN2 Rabbit Monoclonal Antibody**Catalog #: AMRe84595**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200
Molecular Weight	Calculated MW: 9 kDa ; Observed MW: 18 kDa

Antigen Information


Gene Name	HMGN2
Alternative Names	HMG17; HMGN2;;HMGN2
Gene ID	
SwissProt ID	P05204
Immunogen	A synthesized peptide derived from human HMGN2

Background

Binds to the inner side of the nucleosomal DNA thus altering the interaction between the DNA and the histone octamer. May be involved in the process which maintains transcribable genes in a unique chromatin conformation.

Research Area

Image Data



Western blot analysis of HMGN2 expression in K562 cell lysate.