

Product Name: EEF1G Rabbit Monoclonal Antibody**Catalog #: AMRe84599**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,ICC,FC,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.71mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
Molecular Weight	Calculated MW: 50 kDa ; Observed MW: 48 kDa

Antigen Information


Gene Name	EEF1G
Alternative Names	eEF 1B gamma; EEF1G; EF 1G; EF1 gamma; EF1G; GIG35;;eEF1B gamma
Gene ID	
SwissProt ID	P26641
Immunogen	A synthesized peptide derived from human eEF1B gamma

Background


Probably plays a role in anchoring the complex to other cellular components.

Research Area

Image Data



Western blot analysis of EEF1G expression in (1) 293T cell lysate; (2) HepG2 cell lysate; (3) NIH/3T3 cell lysate.



Immunoprecipitate (IP) analysis using the Antibody at 1:50 dilution. (wb at 1:3K dilution)