

Product Name: Glycerol kinase Rabbit Monoclonal Antibody**Catalog #: AMRe84625**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,FC,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,FC 1:20-1:100,IP 1:20-1:50
Molecular Weight	Calculated MW: 61 kDa ; Observed MW: 50 kDa

Antigen Information


Gene Name	Glycerol kinase
Alternative Names	GK1; GKD; Glycerokinase;;Glycerol kinase
Gene ID	
SwissProt ID	P32189
Immunogen	A synthesized peptide derived from human Glycerol kinase

Background

Key enzyme in the regulation of glycerol uptake and metabolism.

Research Area

Image Data



Western blot analysis of Glycerol kinase expression in HepG2 cell lysate.