
Product Name: CTNNA3 Rabbit Monoclonal Antibody**Catalog #: AMRe84628**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200,IP 1:20-1:50
Molecular Weight	100 kDa

Antigen Information


Gene Name	CTNNA3
Alternative Names	Alpha T catenin; Cadherin associated protein; Catenin alpha 3; Catna3; CTNNA3; VR22;;Catenin alpha 3
Gene ID	
SwissProt ID	Q9UI47
Immunogen	A synthesized peptide derived from human Catenin alpha 3

Background

May be involved in formation of stretch-resistant cell-cell adhesion complexes.

Research Area

Image Data



Western blot analysis of CTNNA3 expression in 293 cell lysate.