
Product Name: Histone H3(formylK122) Rabbit Monoclonal Antibody**Catalog #: AMRe84637**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200
Molecular Weight	15 kDa

Antigen Information

Gene Name	Histone H3(formylK122)
Alternative Names	Histone H3;;Formyl-Histone H3 (K123)
Gene ID	
SwissProt ID	P68431
Immunogen	A synthesized peptide derived from human Histone H3.1 around the formylation site of K123


Background

Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA

repair, DNA replication and chromosomal stability. DNA accessibility is regulated via a complex set of post-translational modifications of histones, also called histone code, and nucleosome remodeling.

Research Area

Image Data



Western blot analysis of Histone H3 (formyl K122) expression in HeLa cell lysate.