
Product Name: HLA-DRB4 Rabbit Monoclonal Antibody**Catalog #: AMRe84643**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200,FC 1:20-1:100
Molecular Weight	30 kDa

Antigen Information


Gene Name	HLA-DRB4
Alternative Names	DR4; DRB4; HLA class II histocompatibility antigen; HLA-DRB4;;HLA DRB4
Gene ID	
SwissProt ID	P13762
Immunogen	A synthesized peptide derived from human HLA DRB4

Background

Binds peptides derived from antigens that access the endocytic route of antigen presenting cells (APC) and presents them on the cell surface for recognition by the CD4 T-cells.

Research Area

Image Data



Western blot analysis of HLA-DRB4 expression in Raji cell lysate.