

Product Name: Desmocollin 3 Rabbit Monoclonal Antibody
Catalog #: AMRe84644

Summary

Production Name	Desmocollin 3 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,ICC,FC
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	Desmocollin 3
Alternative Names	CDHF3; Desmocollin3; Desmocollin4; DSC; DSC1; DSC2; DSC3; DSC4; HTCP;;Desmocollin 3
Gene ID	
SwissProt ID	Q14574. A synthesized peptide derived from human Desmocollin 3

Application

Dilution Ratio	WB:1:1000-1:2000,ICC:1:50-1:200,FC:1:20-1:100
Molecular Weight	100 kDa

Product Name: Desmocollin 3 Rabbit Monoclonal Antibody
Catalog #: AMRe84644




Background

Component of intercellular desmosome junctions. Involved in the interaction of plaque proteins and intermediate filaments mediating cell-cell adhesion. May contribute to epidermal cell positioning (stratification) by mediating differential adhesiveness between cells that express different isoforms.

Research Area

Image Data



Western blot analysis of Desmocollin 3 expression in A431 cell lysate.

Note

For research use only.