

Product Name: MRPL28 Rabbit Monoclonal Antibody
Catalog #: AMRe84645



Summary

Production Name	MRPL28 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP,ICC,FC
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	MRPL28
Alternative Names	L28mt; MAAT1; mrpl28; P15;;MRPL28
Gene ID	
SwissProt ID	Q13084. A synthesized peptide derived from human MRPL28

Application

Dilution Ratio	WB:1:1000-1:2000,IP:1:20-1:50,ICC:1:50-1:200,FC:1:20-1:100
Molecular Weight	30 kDa

Background

Potentially represents an important therapeutic reagent for HLA-A24 patients. This antigen is recognized by tumor-


Product Name: MRPL28 Rabbit Monoclonal Antibody
Catalog #: AMRe84645



infiltrating lymphocyte (TIL) 1290 in the context of HLA-A24.

Research Area

Image Data



Western blot analysis of MRPL28 expression in 293 cell lysate.

Note

For research use only.