
Product Name: UCHL3 Rabbit Monoclonal Antibody**Catalog #: AMRe84647**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200
Molecular Weight	26 kDa

Antigen Information


Gene Name	UCHL3
Alternative Names	Ubiquitin thioesterase L3; UCHL3;;UCHL3
Gene ID	
SwissProt ID	P15374
Immunogen	A synthesized peptide derived from human UCHL3

Background

Deubiquitinating enzyme (DUB) that controls levels of cellular ubiquitin through processing of ubiquitin precursors and ubiquitinated proteins.

Research Area

Image Data



Western blot analysis of UCHL3 expression in K562 cell lysate.