

Product Name: AP2S1 Rabbit Monoclonal Antibody**Catalog #: AMRe84651**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.71mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200,IP 1:20-1:50
Molecular Weight	17 kDa

Antigen Information


Gene Name	AP2S1
Alternative Names	AP17; Ap2s1; CLAPS2; Sigma adaptin 3b; Sigma2 adaptin;;AP2S1
Gene ID	
SwissProt ID	P53680
Immunogen	A synthesized peptide derived from human AP2S1

Background

Component of the adaptor protein complex 2 (AP-2). Adaptor protein complexes function in protein Transport via Transport vesicles in different membrane traffic pathways.

Research Area

Image Data



Western blot analysis of AP2S1 expression in (1) 293 cell lysate; (2) Mouse brain lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

