

Product Name: Olig3 Rabbit Monoclonal Antibody**Catalog #: AMRe84659**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000, ICC 1:50-1:200
Molecular Weight	29 kDa

Antigen Information


Gene Name	Olig3
Alternative Names	bHLHB7; bHLHe20; Olig3; Oligo3;;Oligo3
Gene ID	
SwissProt ID	Q7RTU3
Immunogen	A synthesized peptide derived from human Oligo3

Background

May determine the distinct specification program of class A neurons in the dorsal part of the spinal cord and suppress specification of class B neurons.

Research Area

Image Data



Western blot analysis of Olig3 expression in 293T cell lysate.