

Product Name: USP24 Rabbit Monoclonal Antibody**Catalog #: AMRe84661**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	294 kDa

Antigen Information

Gene Name	USP24
Alternative Names	Usp24;;USP24
Gene ID	
SwissProt ID	Q9UPU5
Immunogen	A synthesized peptide derived from human USP24


Background

Ubiquitin-specific protease that regulates cell survival in various contexts through modulating the protein stability of some of its substrates including DDB2, MCL1 or TP53. Plays a positive role on ferritinophagy where ferritin is degraded in lysosomes and

releases free iron.

Research Area

Image Data



Western blot analysis of USP24 expression in HeLa cell lysate.