
Product Name: SNRPA1 Rabbit Monoclonal Antibody**Catalog #: AMRe84665**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200
Molecular Weight	28 kDa

Antigen Information


Gene Name	SNRPA1
Alternative Names	Lea1; SNRPA1;;SNRPA1
Gene ID	
SwissProt ID	P09661
Immunogen	A synthesized peptide derived from human SNRPA1

Background

Involved in pre-mRNA splicing as component of the spliceosome. Associated with sn-RNP U2, where it contributes to the binding of stem loop IV of U2 snRNA.

Research Area

Image Data



Western blot analysis of SNRPA1 expression in HeLa cell lysate.