
Product Name: SUOX Rabbit Monoclonal Antibody**Catalog #: AMRe84668**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.71mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200
Molecular Weight	60 kDa

Antigen Information


Gene Name	SUOX
Alternative Names	Sulfite oxidase; Suox;;SUOX
Gene ID	
SwissProt ID	P51687
Immunogen	A synthesized peptide derived from human SUOX

Background


Catalyzes the oxidation of sulfite to sulfate, the terminal reaction in the oxidative degradation of sulfur-containing amino acids.

Research Area

Image Data



Western blot analysis of SUOX expression in 293T cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.