

Product Name: TUB1 Rabbit Monoclonal Antibody
Catalog #: AMRe84682



Summary

Production Name	TUB1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	TUB1
Alternative Names	rd5; TUB 1; TUB;;Tubby protein homolog
Gene ID	
SwissProt ID	P50607. A synthesized peptide derived from human Tubby protein homolog

Application

Dilution Ratio	WB:1:1000-1:2000,IHC:1:100-1:200
Molecular Weight	Calculated MW: 56 kDa ; Observed MW: 60 kDa

Background

Functions in signal transduction from heterotrimeric G protein-coupled receptors. Binds to membranes containing


Product Name: TUB1 Rabbit Monoclonal Antibody
Catalog #: AMRe84682



phosphatidylinositol 4,5-bisphosphate.

Research Area

Image Data



Western blot analysis of TUB1 expression in Neuro-2a cell lysate.

Note

For research use only.