

Product Name: ACAT2 Rabbit Monoclonal Antibody**Catalog #: AMRe84686**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200
Molecular Weight	Calculated MW: 41 kDa ; Observed MW: 38 kDa

Antigen Information


Gene Name	ACAT2
Alternative Names	Acat2; ACTL;;ACAT2
Gene ID	
SwissProt ID	Q9BWD1
Immunogen	A synthesized peptide derived from human ACAT2

Background


Involved in the biosynthetic pathway of cholesterol.

Research Area

Image Data



Western blot analysis of ACAT2 expression in K562 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.