

Product Name: SRP19 Rabbit Monoclonal Antibody**Catalog #: AMRe84693**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,ICC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
Molecular Weight	Calculated MW: 16 kDa ; Observed MW: 20 kDa

Antigen Information

Gene Name	SRP19
Alternative Names	SRP19;;SRP19
Gene ID	
SwissProt ID	P09132
Immunogen	A synthesized peptide derived from human SRP19


Background

Component of the signal recognition particle (SRP) complex, a ribonucleoprotein complex that mediates the cotranslational targeting of secretory and membrane proteins to the endoplasmic reticulum (ER). Binds directly to 7SL RNA. Mediates binding

of SRP54 to the SRP complex.

Research Area

Image Data



Western blot analysis of SRP19 expression in HeLa cell lysate.