

Product Name: RASA1 Rabbit Monoclonal Antibody**Catalog #: AMRe84695**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200
Molecular Weight	Calculated MW: 116 kDa ; Observed MW: 140 kDa

Antigen Information


Gene Name	RASA1
Alternative Names	CMAVM; GAP; p120GAP; p120RASGAP; PKWS; RASA; RASA1; RasGAP;;RASA1
Gene ID	
SwissProt ID	P20936
Immunogen	A synthesized peptide derived from human RASA1

Background


Inhibitory regulator of the Ras-cyclic AMP pathway. Stimulates the GTPase of normal but not oncogenic Ras p21.

Research Area

Image Data



Western blot analysis of RASA1 expression in Rat brain cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.