
Product Name: PPP6C Rabbit Monoclonal Antibody**Catalog #: AMRe84703**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,ICC 1:50-1:200,FC 1:20-1:100
Molecular Weight	Calculated MW: 35 kDa ; Observed MW: 32 kDa

Antigen Information


Gene Name	PPP6C
Alternative Names	PP6; PP6C; PPP6; Ppp6c;;PPP6C
Gene ID	
SwissProt ID	O00743
Immunogen	A synthesized peptide derived from human PPP6C

Background

Catalytic subunit of protein phosphatase 6 (PP6).PP6 is a component of a signaling pathway regulating cell cycle progression in response to IL2 receptor stimulation.

Research Area

Image Data



Western blot analysis of PPP6C expression in HeLa cell lysate.