

Product Name: DPP3 Rabbit Monoclonal Antibody**Catalog #: AMRe84705**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000, FC 1:20-1:100
Molecular Weight	83 kDa

Antigen Information

Gene Name	DPP3
Alternative Names	Dipeptidyl arylamidase III; Dipeptidyl peptidase 3; Dipeptidyl peptidase III; dpp3; DPPIII;;DPP3
Gene ID	
SwissProt ID	Q9NY33
Immunogen	A synthesized peptide derived from human DPP3


Background

Cleaves and degrades bioactive peptides, including angiotensin, Leu-enkephalin and Met-enkephalin. Also cleaves Arg-Arg-

beta-naphthylamide.

Research Area

Image Data



Western blot analysis of DPP3 expression in 293T cell lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

