

**Product Name: Gamma glutamyl hydrolase Rabbit
Monoclonal Antibody
Catalog #: AMRe84712**

Summary

Production Name	Gamma glutamyl hydrolase Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	Gamma glutamyl hydrolase
Alternative Names	folylpolygamma glutamyl hydrolase; gamma glutamyl hydrolase; GGH; ;gamma Glutamyl hydrolase
Gene ID	
SwissProt ID	Q92820.A synthesized peptide derived from human gamma Glutamyl hydrolase

Application

Dilution Ratio	WB:1:1000-1:2000
Molecular Weight	36 kDa

**Product Name: Gamma glutamyl hydrolase Rabbit
Monoclonal Antibody
Catalog #: AMRe84712**




Background


Hydrolyzes the polyglutamate sidechains of pteroylpolyglutamates. Progressively removes gamma-glutamyl residues from pteroylpoly-gamma-glutamate to yield pteroyl-alpha-glutamate (folic acid) and free glutamate. May play an important role in the bioavailability of dietary pteroylpolyglutamates and in the metabolism of pteroylpolyglutamates and antifolates.

Research Area

Image Data




All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

**Product Name: Gamma glutamyl hydrolase Rabbit
Monoclonal Antibody
Catalog #: AMRe84712**



Western blot analysis of Gamma glutamyl hydrolase expression in MCF7 cell lysate.

Note

For research use only.