
Product Name: MRPS15 Rabbit Monoclonal Antibody**Catalog #: AMRe84719**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,ICC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
Molecular Weight	Calculated MW: 30 kDa ; Observed MW: 22 kDa

Antigen Information


Gene Name	MRPS15
Alternative Names	DC37; MRPS15; RPMS15; S15mt;;MRPS15
Gene ID	
SwissProt ID	P82914
Immunogen	A synthesized peptide derived from human MRPS15

Background


Belongs to the ribosomal protein S15P family.

Research Area

Image Data



Western blot analysis of MRPS15 expression in HeLa cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.