

Product Name: Thromboxane A2 receptor Rabbit Monoclonal Antibody**Catalog #: AMRe84722**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	Calculated MW: 37 kDa ; Observed MW: 50 kDa

Antigen Information


Gene Name	Thromboxane A2 receptor
Alternative Names	BDPLT13; Tbx2r; TP receptor; TXA2R;;TBXA2R
Gene ID	
SwissProt ID	P21731
Immunogen	A synthesized peptide derived from human TBXA2R

Background

Receptor for thromboxane A2 (TXA2), a potent stimulator of platelet aggregation.

Research Area

Image Data



Western blot analysis of Thromboxane A2 receptor expression in human platelet cell lysate.