

Product Name: Syntaxin16 Rabbit Monoclonal Antibody**Catalog #: AMRe84724**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,ICC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
Molecular Weight	37 kDa

Antigen Information


Gene Name	Syntaxin16
Alternative Names	hsyn16; Stx16; Syn16;;Syntaxin 16
Gene ID	
SwissProt ID	O14662
Immunogen	A synthesized peptide derived from human Syntaxin 16

Background

SNARE involved in vesicular transport from the late endosomes to the trans-Golgi network.

Research Area

Image Data



Western blot analysis of Syntaxin 16 expression in HeLa cell lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

