
Product Name: TOE1 Rabbit Monoclonal Antibody**Catalog #: AMRe84731**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	Calculated MW: 57 kDa ; Observed MW: 60 kDa

Antigen Information


Gene Name	TOE1
Alternative Names	Toe1;;Target of EGR1 protein 1
Gene ID	
SwissProt ID	Q96GM8
Immunogen	A synthesized peptide derived from human Target of EGR1 protein 1

Background

Inhibits cell growth rate and cell cycle. Induces CDKN1A expression as well as TGF-beta expression. Mediates the inhibitory growth effect of EGR1. Involved in the maturation of snRNAs and snRNA 3'-tail processing.

Research Area

Image Data



Western blot analysis of TOE1 expression in 293T cell lysate.