
Product Name: Calcipressin 1 Rabbit Monoclonal Antibody**Catalog #: AMRe84733**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	28 kDa

Antigen Information


Gene Name	Calcipressin 1
Alternative Names	Adapt78; Calcipressin 1; CSP1; DSC1; DSCR1; MCIP1; RCAN1; RCN1;;Calcipressin 1
Gene ID	
SwissProt ID	P53805
Immunogen	A synthesized peptide derived from human Calcipressin 1

Background


Inhibits calcineurin-dependent transcriptional responses by binding to the catalytic domain of calcineurin A. Could play a role during central nervous system development.

Research Area

Image Data



Western blot analysis of Calcipressin 1 expression in mouse heart cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.