

---

**Product Name: Sterol carrier protein 2 Rabbit Monoclonal Antibody****Catalog #: AMRe84735**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC,ICC/IF,ICC
<b>Reactivity</b>	Human,Mouse,Rat
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
<b>Molecular Weight</b>	Calculated MW: 59 kDa ; Observed MW: 58 kDa

**Antigen Information**


<b>Gene Name</b>	Sterol carrier protein 2
<b>Alternative Names</b>	NLTP; SCP2; SCPchi; SCPX;;Sterol carrier protein 2
<b>Gene ID</b>	
<b>SwissProt ID</b>	P22307
<b>Immunogen</b>	A synthesized peptide derived from human Sterol carrier protein 2

**Background**

Mediates in vitro the transfer of all common phospholipids, cholesterol and gangliosides between membranes. May play a role in regulating steroidogenesis.

## Research Area

## Image Data



Western blot analysis of Sterol carrier protein 2 expression in HepG2 cell lysate.