

Product Name: UFM1 Rabbit Monoclonal Antibody
Catalog #: AMRe84740



Summary

Production Name	UFM1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP,IHC
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	UFM1
Alternative Names	BM002; Chromosome 13 open reading frame 20; Ubiquitin-fold modifier 1; UFM1;;UFM1
Gene ID	
SwissProt ID	P61960.A synthesized peptide derived from human UFM1

Application

Dilution Ratio	WB:1:1000-1:2000,IP:1:20-1:50,IHC:1:100-1:200
Molecular Weight	Calculated MW: 9 kDa ; Observed MW: 10 kDa

Product Name: UFM1 Rabbit Monoclonal Antibody
Catalog #: AMRe84740




Background


Ubiquitin-like modifier protein which binds to a number of target proteins, such as DDRGK1.

Research Area


Image Data



Western blot analysis of UFM1 expression in HepG2 cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Immunoprecipitate (IP) analysis using the Antibody at 1:50 dilution. (wb at 1:500 dilution)

Product Name: UFM1 Rabbit Monoclonal Antibody
Catalog #: AMRe84740



Note

For research use only.