

Product Name: ARFRP1/ARP Rabbit Monoclonal Antibody**Catalog #: AMRe84751**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,ICC 1:50-1:200
Molecular Weight	23 kDa

Antigen Information


Gene Name	ARFRP1/ARP
Alternative Names	ARFRP1; ARL18; ARP; ARP1;;ARP
Gene ID	
SwissProt ID	Q13795
Immunogen	A synthesized peptide derived from human ARP

Background

Possibly involved in plasma membrane-related signaling events.

Research Area

Image Data



Western blot analysis of ARFRP1 / ARP expression in 293 cell lysate.