
Product Name: Cystatin B Rabbit Monoclonal Antibody**Catalog #: AMRe84753**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200
Molecular Weight	Calculated MW: 11 kDa ; Observed MW: 13 kDa

Antigen Information


Gene Name	Cystatin B
Alternative Names	CHROW21; CST6; CSTB; Cystatin B (stefin B); CYTB EPM1; EPM1A; PME; STFB; ULD;;Cystatin B
Gene ID	
SwissProt ID	P04080
Immunogen	A synthesized peptide derived from human Cystatin B

Background

This is an intracellular thiol proteinase inhibitor. Tightly binding reversible inhibitor of cathepsins L, H and B.

Research Area

Image Data



Western blot analysis of Cystatin B expression in U87-MG cell lysate.