
Product Name: HDLBP Rabbit Monoclonal Antibody**Catalog #: AMRe84762**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200
Molecular Weight	Calculated MW: 141 kDa ; Observed MW: 150 kDa

Antigen Information


Gene Name	HDLBP
Alternative Names	HBP; HDL binding protein; HDLBP; PRO2900; VGL; Vigilin;;Vigilin
Gene ID	
SwissProt ID	Q00341
Immunogen	A synthesized peptide derived from human Vigilin

Background

Appears to play a role in cell sterol metabolism. It may function to protect cells from over-accumulation of cholesterol.

Research Area

Image Data



Western blot analysis of HDLBP expression in HeLa cell lysate.