

Product Name: RASSF2 Rabbit Monoclonal Antibody
Catalog #: AMRe84771



Summary

Production Name	RASSF2 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	RASSF2
Alternative Names	CENP 34; RASFADIN; RASSF2;;RASSF2
Gene ID	
SwissProt ID	P50749. A synthesized peptide derived from human RASSF2

Application

Dilution Ratio	WB:1:1000-1:2000
Molecular Weight	Calculated MW: 38 kDa ; Observed MW: 37 kDa

Background

Potential tumor suppressor. Acts as a KRAS-specific effector protein. May promote apoptosis and cell cycle arrest. Stabilizes


Product Name: RASSF2 Rabbit Monoclonal Antibody
Catalog #: AMRe84771



STK3/MST2 by protecting it from proteasomal degradation.

Research Area

Image Data



Western blot analysis of RASSF2 expression in Jurkat cell lysate.

Note

For research use only.