

Product Name: MVK Rabbit Monoclonal Antibody**Catalog #: AMRe84772**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000, FC 1:20-1:100
Molecular Weight	42 kDa

Antigen Information


Gene Name	MVK
Alternative Names	LRBP; Mevalonate kinase 1; mvk; MVLK; POROK3;; Mevalonate kinase
Gene ID	
SwissProt ID	Q03426
Immunogen	A synthesized peptide derived from human Mevalonate kinase

Background

Catalyzes the phosphorylation of mevalonate to mevalonate 5-phosphate, a key step in isoprenoid and cholesterol biosynthesis.

Research Area

Image Data



Western blot analysis of MVK expression in HepG2 cell lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

