
Product Name: RBM3 Rabbit Monoclonal Antibody**Catalog #: AMRe84775**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,ICC,FC
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100
Molecular Weight	17 kDa

Antigen Information

Gene Name	RBM3
Alternative Names	Rbm3; RNPL;;RBM3
Gene ID	
SwissProt ID	P98179
Immunogen	A synthesized peptide derived from human RBM3


Background

Cold-inducible mRNA binding protein that enhances global protein synthesis at both physiological and mild hypothermic temperatures. Reduces the relative abundance of microRNAs, when overexpressed. Enhances phosphorylation of translation

initiation factors and active polysome formation.

Research Area

Image Data



Western blot analysis of RBM3 expression in Jurkat cell lysate.