

Product Name: PLSCR3 Rabbit Monoclonal Antibody
Catalog #: AMRe84779



Summary

Production Name	PLSCR3 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	PLSCR3
Alternative Names	PLS3; Plscr3;;PLSCR3
Gene ID	
SwissProt ID	Q9NRY6. A synthesized peptide derived from human PLSCR3

Application

Dilution Ratio	WB:1:1000-1:2000
Molecular Weight	32 kDa

Background

Catalyzes calcium-induced ATP-independent rapid bidirectional and non-specific movement of the phospholipids (lipid


Product Name: PLSCR3 Rabbit Monoclonal Antibody
Catalog #: AMRe84779



scrambling or lipid flip-flop) between the inner and outer membrane of the mitochondria. Plays an important role in mitochondrial respiratory function, morphology, and apoptotic response.

Research Area

Image Data



Western blot analysis of PLSCR3 expression in BxPC 3 cell lysate.

Note

For research use only.