
Product Name: EXOC3 Rabbit Monoclonal Antibody**Catalog #: AMRe84791**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,ICC,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,IP 1:20-1:50
Molecular Weight	Calculated MW: 86 kDa ; Observed MW: 80 kDa

Antigen Information


Gene Name	EXOC3
Alternative Names	Exoc3; Sec6; SEC6L1; Sec6p;;Exocyst complex component 3
Gene ID	
SwissProt ID	O60645
Immunogen	A synthesized peptide derived from human Exocyst complex component 3

Background

Component of the exocyst complex involved in the docking of exocytic vesicles with fusion sites on the plasma membrane.

Research Area

Image Data



Western blot analysis of EXOC3 expression in MCF7 cell lysate.