
Product Name: PNK/PNKP Rabbit Monoclonal Antibody**Catalog #: AMRe84797**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,ICC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
Molecular Weight	57 kDa

Antigen Information


Gene Name	PNK/PNKP
Alternative Names	DEM1; EIEE10; MCSZ; PNK1; Pnkp;;PNKP
Gene ID	
SwissProt ID	Q96T60
Immunogen	A synthesized peptide derived from human PNKP

Background

Plays a key role in the repair of DNA damage, functioning as part of both the non-homologous end-joining (NHEJ) and base excision repair (BER) pathways.

Research Area

Image Data



Western blot analysis of PNK / PNKP expression in Jurkat cell lysate.