

Product Name: TSFM Rabbit Monoclonal Antibody**Catalog #: AMRe84803**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,ICC 1:50-1:200,IP 1:20-1:50
Molecular Weight	Calculated MW: 35 kDa ; Observed MW: 32 kDa

Antigen Information


Gene Name	TSFM
Alternative Names	COXPD3; EFTS; EFTSMT; TSFM; TSMT;;TSFM
Gene ID	
SwissProt ID	P43897
Immunogen	A synthesized peptide derived from human TSFM

Background

Associates with the EF-Tu.GDP complex and induces the exchange of GDP to GTP. It remains bound to the aminoacyl-tRNA.EF-Tu.GTP complex up to the GTP hydrolysis stage on the ribosome.

Research Area

Image Data



Western blot analysis of TSFM expression in Jurkat cell lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

