

**Product Name: Phospho-Glucocorticoid Receptor
(Ser226) Rabbit Monoclonal Antibody**
Catalog #: AMRe84904

Summary

Production Name	Phospho-Glucocorticoid Receptor (Ser226) Rabbit Monoclonal Antibody
Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Phosphorylated
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	Phospho-Glucocorticoid Receptor (Ser226)
Alternative Names	NR3C1; GRL; Glucocorticoid receptor; GR; Nuclear receptor subfamily 3 group C member 1
Gene ID	2908.0
SwissProt ID	P04150.A synthetic phosphopeptide corresponding to residues surrounding Ser226 of human Glucocorticoid Receptor

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100
Molecular Weight	Calculated MW: 86 kDa; Observed MW: 94 kDa

**Product Name: Phospho-Glucocorticoid Receptor
(Ser226) Rabbit Monoclonal Antibody**
Catalog #: AMRe84904




Background

The protein encoded by this gene is a receptor for glucocorticoids and can act as both a transcription factor and a regulator of other transcription factors. The encoded protein can bind DNA as a homodimer or as a heterodimer with another protein such as the retinoid X receptor. This protein can also be found in heteromeric cytoplasmic complexes along with heat shock factors and immunophilins.


Research Area

Apoptosis

Image Data



Western blot analysis of Phospho-Glucocorticoid Receptor (Ser226) in HeLa, A549, C6 lysates using Phospho-Glucocorticoid Receptor (Ser226) antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Glucocorticoid Receptor (Phospho-Ser226) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



**Product Name: Phospho-Glucocorticoid Receptor
(Ser226) Rabbit Monoclonal Antibody
Catalog #: AMRe84904**

Note

For research use only.