
Product Name: Phospho-TAOK1/2/3 (Ser181/Ser181/Ser177) Rabbit Monoclonal Antibody
Catalog #: AMRe84919

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Phosphorylated
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100
Molecular Weight	Calculated MW: 116 kDa; Observed MW: 105,116,138 kDa

Antigen Information

Gene Name	Phospho-TAOK1/2/3 (Ser181/Ser181/Ser177)
Alternative Names	DPK; JIK; hKFC-A; MAP3K18
Gene ID	51347/57551/9344
SwissProt ID	Q7L7X3/Q9H2K8/Q9UL54
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Ser181 of human TAOK1


Background

Serine/threonine-protein kinase involved in various processes such as p38/MAPK14 stress-activated MAPK cascade, DNA damage response and regulation of cytoskeleton stability. Phosphorylates MAP2K3, MAP2K6 and MARK2. Acts as an activator


of the p38/MAPK14 stress-activated MAPK cascade by mediating phosphorylation and subsequent activation of the upstream MAP2K3 and MAP2K6 kinases. Involved in G-protein coupled receptor signaling to p38/MAPK14. In response to DNA damage, involved in the G2/M transition DNA damage checkpoint by activating the p38/MAPK14 stress-activated MAPK cascade, probably by mediating phosphorylation of MAP2K3 and MAP2K6.

Research Area


Image Data



Western blot analysis of Phospho-TAOK1(Ser181)/TAOK2 (Ser181)/TAOK3 (Ser177) in 3T3 lysates using Phospho-TAOK1/2/3 (Ser181/Ser181/Ser177) antibody.



Western blot analysis of TAOK1 (Phospho-Ser181)/TAOK2 (Phospho-Ser181)/TAOK3 (Phospho-Ser177) in HeLa, A549, HL-60, U2OS, C6 lysates using TAOK1 (Phospho-Ser181)/TAOK2 (Phospho-Ser181)/TAOK3 (Phospho-Ser177) antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using TAOK1 (Phospho-Ser181)/TAOK2 (Phospho-Ser181)/TAOK3 (Phospho-Ser177) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.