

Product Name: ZNF364 Rabbit Monoclonal Antibody**Catalog #: AMRe85110**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 34 kDa; Observed MW: 34 kDa

Antigen Information


Gene Name	ZNF364
Alternative Names	E3 ubiquitin-protein ligase RNF115; Rabring 7; Zinc finger protein 364
Gene ID	27246.0
SwissProt ID	Q9Y4L5
Immunogen	A synthetic peptide of human ZNF364

Background


E3 ubiquitin-protein ligase that mediates E2-dependent, 'Lys-48'- and/or 'Lys-63'-linked polyubiquitination of substrates and may play a role in diverse biological processes.

Research Area

Image Data



Western blot analysis of ZNF364/RFN115 in PC-12 lysates using ZNF364 antibody.



Western blot analysis of ZNF364/RFN115 in Jurkat, HeLa lysates using ZNF364/RFN115 antibody.