
Product Name: TAF3 Rabbit Monoclonal Antibody**Catalog #: AMRe85124**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100
Molecular Weight	Calculated MW: 104 kDa; Observed MW: 140 kDa

Antigen Information


Gene Name	TAF3
Alternative Names	Transcription initiation factor TFIID subunit 3
Gene ID	83860.0
SwissProt ID	Q5VWG9
Immunogen	Recombinant protein of human TAF3

Background


The TFIID basal transcription factor complex plays a major role in the initiation of RNA polymerase II (Pol II)-dependent transcription.

Research Area

Image Data



Western blot analysis of TAF3 in Jurkat lysates using TAF3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using TAF3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.