

Product Name: ARFRP1 Rabbit Monoclonal Antibody**Catalog #: AMRe85130**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000, IP 1:10-1:20
Molecular Weight	Calculated MW: 23 kDa; Observed MW: 23 kDa

Antigen Information


Gene Name	ARFRP1
Alternative Names	ARP; Arp1; ARL18
Gene ID	10139.0
SwissProt ID	Q13795
Immunogen	A synthetic peptide of human ARFRP1

Background

Trans-Golgi-associated GTPase that regulates protein sorting.

Research Area

Image Data



Western blot analysis of ARFRP1/ARP in C6 lysates using ARFRP1 antibody.