
Product Name: OXGR1 Rabbit Monoclonal Antibody**Catalog #: AMRe85132**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 38 kDa; Observed MW: 38 kDa

Antigen Information


Gene Name	OXGR1
Alternative Names	P2Y15; Alpha-ketoglutarate receptor 1; G-protein coupled receptor 80; G-protein coupled receptor 99; P2Y purinoceptor 15
Gene ID	27199.0
SwissProt ID	Q96P68
Immunogen	A synthetic peptide of human OXGR1

Background


Receptor for alpha-ketoglutarate.

Research Area

Image Data



Western blot analysis of OXGR1/GPR99 in rat colon lysates using OXGR1 antibody.



Western blot analysis of OXGR1/GPR99 in Hela lysates using OXGR1/GPR99 antibody.