

Product Name: CENPA Rabbit Monoclonal Antibody**Catalog #: AMRe85136**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 16 kDa; Observed MW: 16 kDa

Antigen Information


Gene Name	CENPA
Alternative Names	Centromere autoantigen A; Centromere protein A; CENP-A
Gene ID	1058.0
SwissProt ID	P49450
Immunogen	A synthetic peptide of human CENPA

Background

Histone H3-like nucleosomal protein that is specifically found in centromeric nucleosomes.

Research Area

Image Data



Western blot analysis of CENPA in K562, HeLa lysates using CENPA antibody.